

PR2 Maintenance

ROS + PR2 Training Workshop


PR2 Maintenance Overview

- PR2 overview
- Review electrical and ESD safety
- Review PR2 power system
- PR2 accessory kit
- PR2 tool kit
- Hands on demos
- See *support.willowgarage.com* to order replacement parts


PR2 Overview

- PR2 “anatomy”
- Hands-on demos
 - Filter maintenance
 - Gripper tip replacement
 - Counter balance adjustment
- Points of interest
 - Power switch
 - A/C breaker
 - Run-Stop button
 - Plug for the ESD strap


PR2 Front View



PR2 Side View



PR2 Rear View



Electrical and ESD safety

- When performing maintenance:
 - Turn off power even when removing covers
 - The only exceptions are when performing counter balance adjustments and replacing the fan filter
 - Disconnect AC power to the PR2
 - Always use an ESD strap
- ESD straps are provided in the accessory kit:
 - Place strap around bare wrist
 - Plug jack into accessory panel

PR2 Power System

PR2 internal power precautions and safety features:

- Battery power is “live” within the robot, *even when the power is off*
- The power board controls all the power to the PR2 systems and is accessible/controllable via the ROS `pr2_power_node`
- There are as many as 7 levels of protection but you cannot assume they will prevent damage or injury
- Always turn off the PR2 power when performing maintenance except when specifically required to complete a procedure such as the counter balance adjustment

PR2 Power Requirements

Wall power:

- The PR2 can only be plugged into a dedicated 15A ***minimum*** circuit
- Only use the supplied power cord
- Only use the supplied self plug-in cord
- See support.willowgarage.com for ordering replacement cords

Battery runtime

- Approximately 2 hours stationary, less when moving

PR2 Accessory Kit

15 FT power cord 68-06167

Self plug-in cord 68-06153

Joystick 68-02985

Hard drive bay keys 68-03337

Wireless run-stop 68-01059

USB joystick cable 68-03533

Small Checkerboard 68-03823

Fan filters 68-02802 (x5)

ESD strap 68-03024

Base station recovery DVD

User manual

Gripper tips box:

Sensor fingertips 68-03840
(x5) (use only with
boots!)

Fingertip boots 68-03347
(x20)



PR2 Tool Kit (1)

Counter balance adjustment rod -
10mm hex 68-07034

Counter balance adjustment tube
68-07147

Ratchet head bits:

1.5mm 68-07040
2.5mm 68-07042

Driver bits:

1.5mm 68-07062
2.0mm 68-07063
2.5mm 68-07064
2.5mm ball end 68-07059
3.0mm 68-07065
4.0mm 68-07066
5.0mm 68-07067

Needle nose pliers 68-07112

Diagonal Cutter 68-07146

Zip-ties 68-07150

1/4 sq. male to 1/4 hex male 68-07075

1/4 sq. drive to 1/4 hex female 68-07077


3/8 sq. drive to 1/4 hex female 68-07078

10mm 3/8 drive socket 68-07079

13mm 1/4 drive socket 68-07085

13mm 3/8 drive socket 68-07086

16mm 3/8 drive socket 68-07142



PR2 Tool Kit (2)

7/32 1/4 drive socket 68-07090

3/8 drive short extension, 68-07116

3/8 drive wobble extension, 68-07118

1/4 drive quick release extension,
68-07140

1/4 drive flex joint, 68-07141

1/4 drive ratchet wrench, 68-07138

3/8 drive ratchet wrench, 68-07139

Needle Nose pliers, 68-07112

Mountz wrench and driver tools:

Fixed 0.34Nm wrench, with 1.5mm right
angle Allen, 68-07025


Fixed 0.68Nm wrench, 68-07028

Fixed 2.0Nm wrench, 68-07029

1/4 hex ratchet head, 68-07033 (for use
with 0.68 and 2.0 Nm wrenches)

0.20 - 1.20 Nm torque tool, 68-07014
(small)

1 - 6 Nm Torque tool, 68-07015 (large)



Using Torque Tools

Small and large driver

Special stack-ups



See support.willowgarage.com:

PR2 Toolkit: Setting Torque Driver

Blue Adjustment
Nob - Small
Torque Driver




Red Adjustment
Nob - Large
Torque Driver




Using Loctite

- Use Loctite to retain fasteners and keep them from vibrating loose
- There is no need to clean the old Loctite
- The fresh Loctite will integrate with the old



Only a small "patch" of Loctite is needed.



Old Loctite.

Today's Demos

Video of gripper tip replacement

Live demo of fan filter repair/replacement

Hands-on demo of counterbalance adjustment

Willow Garage Support

HOME ROBOT DOCUMENTATION AND FORUMS MANAGE ACCOUNT RECENT VIEWS NEW

Willow Garage Support

Welcome to Willow Garage Support
Melnee Wise Apr 21 • Announcements (unpin)

PR2 Documentation

- PR2 Manual**
Lost your manual? Looking for updates? Click here for the PR2 manual.
- Repair and Replacement Procedures**
Videos and directions for repairing and replacing PR2 components.
- Modularity Specification**
Want to put a custom gripper or sensor on the PR2? Check out the modularity specification.
- Spare Parts**
Need more sensor tips? Click here for information on requesting new parts.

Mailing Lists & Tickets

- ROS Users Mailing List**
ROS users
Want to ask a question about ROS? Sign up for the mailing list.
- PR2 Users Mailing List**
PR2 users
Want to ask a question about PR2? Sign up for the mailing list.
- PR2 Admin Mailing List**
PR2 admin
Are you a PR2 admin looking for help? Sign up for the mailing list.
- PR2 Help Request**
Have a problem with your PR2? File a ticket so we can help you.

See ***support.willowgarage.com*** for all maintenance procedures

Gripper Tip R/R Video

Turn off power



Wear ESD strap

Always cover sensors with boots

Replace boots as soon as they show signs of damage

Willow Garage will not replace damaged sensor tips

See support.willowgarage.com for replacements



Fan Filter Demo

Clean regularly every 1-4 weeks

Clogged filter will cause fan speed to increase
(noisy robot)

Vacuuming filter is not recommended since this
may tear the filter

Maintenance:

- Remove the filter cover plate

- Wash and ***dry*** filter

- Having the fan on helps when placing the filter

- Re-install the cover plate

Counter Balance Review

Overview: PR2 counter balance (CB) holds up the weight of the arms using a spring and pulley system


Goal: The CB should balance the arm so that almost no torque is required to maintain position of the arm

Performance: A well adjusted counter balance will maximize the performance of the PR2 arm controllers

CB Adjustment Overview

With time, the springs of the CB system can start to sag or give. To maintain performance, you will have to *adjust* the mechanism.

- Secondary spring can be stiffened or loosened
 - Clockwise (CW) turns increases force on arm
- “Arm gimbal shaft” controls level arm length of CB spring force
 - Counterclockwise (CCW) increases force on arm



CB Adjustment Demo

- Demonstrate checking and adjusting the counterbalance on the PR2
 - We will demonstrate up front
 - You will adjust your own counterbalance with help from Willow Garage engineers and technicians

Notes:

Safety first

You need two people to adjust the mechanism

Observe ESD protocol (wrist strap)

Watch for pinch points

CB Check Program

- We use a counterbalance adjustment/check program to adjust the counterbalance.
 - The program moves the arm of the robot and measures the torque on each joint
 - It uses training data to recommend adjustments to for a minimum torque setting
 - The program is an iterative solver, in some cases more than one iteration is required

You will run the check program and determine your required adjustment.

CB Check Program Install





Counterbalance adjustment program is *unreleased*.

- You will have to install it from wg-ros-pkg trunk and build it
- Tutorial:
 - ros.org/wiki/qualification/Tutorials
- Follow tutorial instructions to install and build

Important: Do not start the CB adjustment program until instructed to do so

CB Adjustment Demo





Observe safety and ESD protocols.
Remove top covers, Doc6



CB Adjustment Demo

Observe safety and
ESD protocols.

Remove shoulder
covers, Doc16




CB Adjustment Demo

Adjust counter balance, Doc5

Re-run the Check program and adjust the CB if needed,
Doc5

Re-install the covers, Doc6



CB Adjustment Procedure

Check CB with adjustment program

Clear out 2m space around your robot

Turn robot off (*sudo pr2-shutdown*)

Remove covers

Ask your Willow Garage helper for questions

Adjust CB

Make adjustments check tool recommends

Turn robot on

Check CB with adjustment program

Turn robot off

Reinstall covers

Remember!

Safety First

support.willowgarage.com

Always wear an ESD wrist strap

Don't be afraid to ask questions

Questions?????